


Centering: Center LC Media III lenses according to the monocular distance PD and measure at pupil center when looking straight ahead in horizontal (zero) distance direction.



Required ordering data: Refracted prescription for distance and the choice of required support level LC Media III 5, LC Media III 9 or LC Media III 13.

LC Media III are enhanced, everyday single vision lenses; however, they should be blocked and finished similar to progressive lenses.

Materials	50 = 1.50	53 = 1.53	59 = 1.59	60 = 1.60	67 = 1.67
-----------	-----------	-----------	-----------	-----------	-----------

OO = LC Media III	5, 9, 13 = Powerboost	5 = +0.57D	9 = +0.95D	13 = +1.32D
-------------------	-----------------------	------------	------------	-------------